



MARRI LAXMAN REDDY INSTITUTE OF PHARMACY

(Approved by AICTE & PCI, New Delhi and Affiliated to JNTUH)

Dundigal - Gandimaisamma (V) &(M), Medchal (Dt), Hyderabad, Telangana - 500 043.

III/II B. PHARMACY

PRACTICAL MANUAL

OF

MEDICINAL CHEMISTRY-III

ACADEMIC YEAR:2020-2021

NAME: _____

REG.NO: _____

YEAR: _____

About MLRIP



To be an educational Institute of par excellence and produce competent pharmacy professionals to serve the community through research and the ever-increasing needs of Industry.



1. Imparting quality education and innovative research for various career opportunities.
2. Creating conducive academic environment to produce competent pharmacy professionals.
3. Indoctrination of students adorned with high human values and make them aware of their responsibility as health care professionals.

Program Educational Objectives

PEO 1: To produce graduates with sound theoretical knowledge and technical skills required for their career opportunities in various domains.

PEO 2: To incite the students towards research and to address the challenges with their innovative contributions for the benefit of the mankind.

PEO 3: To instill the essence of professionalism, ethical commitment to become a health care professional with sound integrity and adherence to the core human values in the service of the society.



PROGRAM OUTCOMES

1. **Pharmacy Knowledge:** Possess knowledge and comprehension of the core and basic knowledge associated with the profession of pharmacy, including biomedical sciences; pharmaceutical sciences; behavioral, social, and administrative pharmacy sciences; and manufacturing practices.
2. **Planning Abilities:** Demonstrate effective planning abilities including time management, resource management, delegation skills and organizational skills. Develop and implement plans and organize work to meet deadlines.
3. **Problem analysis:** Utilize the principles of scientific enquiry, thinking analytically, clearly and critically, while solving problems and making decisions during daily practice. Find, analyze, evaluate and apply information systematically and shall make defensible decisions.
4. **Modern tool usage:** Learn, select, and apply appropriate methods and procedures, resources, and modern pharmacy-related computing tools with an understanding of the limitations.
5. **Leadership skills:** Understand and consider the human reaction to change, motivation issues, leadership and team-building when planning changes required for fulfillment of practice, professional and societal responsibilities. Assume participatory roles as responsible citizens or leadership roles when appropriate to facilitate improvement in health and well-being.
6. **Professional Identity:** Understand, analyze and communicate the value of their professional roles in society (e.g. health care professionals, promoters of health, educators, managers, employers, employees).
7. **Pharmaceutical Ethics:** Honour personal values and apply ethical principles in professional and social contexts. Demonstrate behavior that recognizes cultural and personal variability in values, communication and lifestyles. Use ethical frameworks; apply ethical principles while making decisions and take responsibility for the outcomes associated with the decisions.
8. **Communication:** Communicate effectively with the pharmacy community and with society at large, such as, being able to comprehend and write effective reports, make effective presentations and documentation, and give and receive clear instructions.
9. **The Pharmacist and society:** Apply reasoning informed by the contextual knowledge to assess societal, health, safety and legal issues and the consequent responsibilities relevant to the professional pharmacy practice.
10. **Environment and sustainability:** Understand the impact of the professional pharmacy solutions in societal and environmental contexts, and demonstrate the knowledge of, and need for sustainable development.
11. **Life-long learning:** Recognize the need for and have the preparation and ability to engage in independent and life-long learning in the broadest context of technological change. Self-assess and use feedback effectively from others to identify learning needs and to satisfy these needs on an ongoing basis.

CONTENT

S. No	Experiment Name	Page No.
1.	Synthesis of Sulphanilamide	01
2.	Synthesis of 7- Hydroxy -4- methyl coumarin	03
3.	Synthesis of Chlorobutanol	05
4.	Synthesis of Tolbutamide	07
5.	Synthesis of Hexamine	09
6.	Assay of Isonicotinic acid hydrazide	11
7.	Assay of Metronidazole	13
8.	Assay of Dapsone	16
9.	Assay of Chlorpheniramine Maleate	18
10.	Assay of Benzyl Penicillin	20
11.	Synthesis of Phenytoin from Benzil by Microwave Irradiation	23
12.	Synthesis of Aspirin Assisted by Microwave Oven	26
13.	Drawing structure and Reaction using Chemdraw	28

Experiment No: 01

Synthesis of Sulphanilamide


Aim:

To synthesis and submit sulphanilamide from p-acetamido benzene sulphanilamide and calculate its percentage yield.

Principle:

Sulphanilamide can be prepared by the reaction of P-acetamido benzene sulphanilamide with Hydrochloric acid or ammonium carbonate. The acetamido groups are easily undergo acid catalysed hydrolysis reaction to form p-amino benzene sulphonamide.

Reaction:



Chemical Required:

Resorcinol	-	1.2 g
Ethyl acetoacetate	-	2.4 ml
Conc. Sulphuric acid	-	7.5 ml

Procedure:

1.5 gm of 4- acetamido benzene sulphonamide is treated with a mixture of 1 ml of conc. Sulphuric acid diluted with 2 ml water. This mixture is gently heated under reflux for 1 hour. Then 3ml of water is added and the solution is boiled again, with the addition of a small quantity of activated charcoal. The solution is filtered while hot, and the filtrate is neutralised with powdered sodium carbonate with stirring until all effervescence ceases and the sulphanilamide is precipitated as a white powder. The solution is cooled,

filtered, the sulphanilamide wash with water and dried. Finally, crude sulphanilamide is recrystallized from hot water, to get colourless crystals.

Melting Point: 163⁰C

Category: Bacteriostatic agent

Report:

The Sulphanilamide was synthesised and submitted, Reported the following

1. Theoretical Yield =
2. Practical Yield =
3. Percentage Yield =

Experiment No: 02

Synthesis of 7- Hydroxy -4- Methyl Coumarin


Aim:

To Synthesis and submit Synthesis of 7- Hydroxy -4- methyl coumarin.

Principle:

A general synthesis of coumarins involves the interaction of a phenol with a β - ketoester in presence of an acid condensing agent (Pechmann reaction). Concentrated sulphuric acid is usually used as the condensing agent for simple monohydric phenols and β - ketoesters, although phenol itself reacts better in the presence of aluminium chloride. The mechanism of the reaction to involve the initial formation of a β - hydroxy ester, which then cyclise and dehydrates to yield the coumarin. Polyhydric phenols, particularly here the two hydroxy groups are meta oriented, react with great ease and sulphuric acid is used as the condensing agent with careful temperature control to ensure a good yield.

Reaction



Chemicals Requirements:

Conc. Sulphuric acid-7.5 ml	Ice – q.s
Resorcinol – 1.2 g	Sodium Hydroxide- q.s
Ethylacetoacetate – 2.4 ml	Ethanol – q.s

Procedure

Take 1 litre of concentrated sulphuric acid in a 3 litres capacity 3 necked flask fitted with a thermometer, mechanical stirrer and a dropping funnel. Immerse the flask in an ice bath. When the temperature falls below 10⁰C, add a solution of 100 g (0.91 mol) of resocinol in 134 g (130.5 ml, 1.03 mol) of redistilled ethyl acetoacetate drop wise and with stirring. Maintain the temperature below 10⁰C by means of an ice-bath during the addition (2h). Keep the mixture at room temperature for 18 hr and then pour it into a mixture of 2 kg of crushed ice with vigours stirring and add 3 litter of water. Collect the

precipitate by suction filtration and wash it with three 25 ml portion of cold water. Dissolve the solid in 1500 ml of 5 percent sodium hydroxide solution, filter and add dilute 2 M sulphuric acid (about 550 ml) with vigourous stirring until the solution is acid to litmus. Filter the crude 4-methyl-7-hydroxy coumarin at the pump, wash with four 25 ml portions of cold water and dry at 100°C. The yield is 155g (97%). Recrystallise from 95 percent ethanol: the pure compound separate in colourless needles, m.p. 185°C.

Use

Laser dye.

Starting material for production of the insecticide ‘hymecromone’.

Report:


7- hydroxy -4- methyl coumarin was synthesized and the percentage yield was found to be%.

Experiment No: 03**Synthesis of Chlorbutanol****Aim:**

To Prepare and submit chlorbutanol from acetone and calculate the Percentage Yield.

Principle:

Chlorbutanol is also known as Chloroketone. It is a trichloro derivative of tertiary butyl alcohol. It is prepared by combination with acetone and chloroform in the presence of solid potassium hydroxide Chlorobutanol is used as local anaesthetic in dental preparation and also antiseptic.

Reaction:**Chemical Required:**

Chloroform – 11 ml

Potassium hydroxide – 1 gm

Acetone – 14 ml

Procedure:

About 11 ml of chloroform and 1 gm of solid potassium hydroxide are taken in a round bottom flask shake until the potassium hydroxide dissolves. To this 14 ml of acetone was added and shaken for 15 minutes. Then set aside for half an hour and crystal of chlorbutanol was separated out.

Melting Point: 95°C – 99°C

Use:

Local anaesthetics

Report:

Chlorbutanol was Prepared and submitted and reported the following

1. Theoretical Yield =
2. Practical Yield =
3. Percentage yield =

Experiment No: 04

Synthesis of Tolbutamide


Aim:

To prepare and submit tolbutamide from p-toluene sulfonamide and calculate its percentage yield.

Principle:

The synthesis of tolbutamide involves addition reaction of p-toluene sulfonamide and butyl isocyanate in the presence of triethylamine and tetrahydrofuran.

Reaction:



Chemical Required:

p-Toluene Sulphonamide – 1 g

Butyl isocyanate – 1 g

Triethylamine – 1.2 ml

Tetrahydrofuran – 10 ml

Procedure:

About N-butyl isocyanate (1m mol) and triethylamine (1.2m mol) in a round bottom flask containing 10 ml of tetrahydrofuran, kept in an ice bath. To the above mixture P-toluene sulfonamide (1m mol) was added drop wise at 0oc. After completing the addition, the temperature was suddenly raised to 35-45oC

and stirred for 3-4 Hrs. Then the solution was filtered. The product was separated and dried. Then it was recrystallised by using Diethyl ether.

Melting Point: 128°C to 129°C

Use:

Anti-diabetic

Report:

Tolbutamide was Prepared and submitted. Report the Following

Theoretical Yield:

Practical Yield:


Percentage Yield:

Experiment No: 05**Synthesis of Hexamine****Aim:**

To prepare and submit hexamine from formaldehyde and calculate its Percentage Yield.

Principle:

Hexamine is heterocyclic organic compound $(\text{CH}_2)_6\text{N}_4$. It has symmetrical tetrahedral cage-like structure. It is prepared by condensation reaction between formaldehyde and ammonia.

Reaction:**Chemical Required:**

Formaldehyde – 4.7 g

Ammonia Solution – 7 g

Procedure:

About 4.7g of 30% formaldehyde solution was taken in a beaker and add 7g of 24% ammonia solution, until the solution is slightly alkaline. The mixture was heated on a water bath for 5 minutes and allowed to stand for 15 minutes. The solution was filtered and then evaporated on a direct flame using china dish to a thick paste. The hexamine crystals are obtained and dried. It was recrystallized from water or alcohol. Hexamine forms colourless, odourless crystals, which are soluble in water and 90% alcohol.

Use:

Urinary anti-infective agent

Report:

Hexamine was prepared and submitted. Report the following

Theoretical Yield:

Practical Yield:

Percentage Yield:

Experiment No: 06

Assay of Isonicotinic acid hydrazide

Aim:



To determine the percentage purity of given sample of Isonicotinic acid hydrazide tablet.

Principle:



Isonicotinic acid hydrazide/ isoniazid is an anti-tubercular drug. It is assayed by the direct titration of Potassium bromate with addition of potassium bromide in the presence of acid medium using hydrochloric acid. During oxidation reaction liberated bromine reacts with isoniazid in aqueous solution to form isonicotinic acid. Azo dye methyl red indicator is used in this titration, decolorized the red colour is the end point.

Reaction:

Standardisation:



Assay:



Procedure:

Preparation of 0.0167 M Potassium bromate solution:

About 5.566gm of Potassium bromate was dissolved in water and made up to 1000 ml with distilled water.

Standardisation of 0.0167 M Potassium bromate:

20 ml of above solution was transferred in a glass stopper flask and 3 g of Potassium iodide and followed with 3 ml of HCl was added. Allow to stand for 5 min, titrate liberated Iodine with 0.1 M sodium thiosulphate adding 3 ml of starch solution TS and the end point is approached.

$$\text{Concentration of 0.0167 M Potassium bromate} = \frac{\text{Molarity of Sodium thiosulphate}}{\text{Volume of Potassium bromate}}$$

Assay of Isonicotinic acid hydrazide tablet

Twenty tablets were weighed accurately and pulverized. A weighed quantity of the tablet power equivalent to 0.25 mg INH was transferred into a clean and dry 100 ml volumetric flask, then sufficient water was added to produce 100 ml. 20 ml of above solution was taken. Then 100 ml of water, 20 ml of hydrochloric acid and 0.2 gm of Potassium bromide were added. Then titrated against slowly with continuous shaking with 0.0167M potassium bromate using 0.05ml of methyl red as an indicator until red colour disappears.

Each ml of 0.0167M Potassium bromate $\text{KBrO}_3 = 0.003439\text{g}$ of $\text{C}_6\text{H}_7\text{N}_3\text{O}$.

Report:

The Molarity of 0.1 M Sodium nitrite

The percentage purity of given isonicotinic acid Hydrazide tablet was found to be:

Experiment No: 07

Assay of Metronidazole

Aim:


To determine the percentage purity of given sample of Metronidazole tablet.

Principle:


Substance with very weak acid or basic property do not give sharp end point in aqueous solution can be titrated with non-aqueous solvents. Non-aqueous titration is based on Bronsted-Lowry theory According to this theory an acid is a proton donor, i.e. a substance which tends to dissociate to yield a proton, and a base is proton acceptor, i.e. a substance which tends to combine with a proton. When an acid HB dissociates it yields a proton together with the conjugate base B of the acid.

Metronidazole is a weak base and it is assayed by non-aqueous titration. When a weakly basic drug is dissolved in acidic solvent with titration of perchloric acid, the basic property of mineral is enhanced. If a strong acid like perchloric acid is added in acetic acid, the acetic acid behaves as a base and combine with proton donated by perchloric acid to form onium ion which act as strong acid, the onium ion is readily donates the proton to the base. The indicator used in this titration is crystal violet. The use of acetic anhydride in this solution, it combines with the water molecule to form acetic acid.


Reaction:



Standardisation:



Assay:



Procedure:

Preparation of Perchloric acid 0.1M

Mixed 8.5 ml of Perchloric acid in 900 ml of glacial acetic acid shake vigorously and continuous stirring. Then added 30 ml acetic anhydride and make up the volume to 1 litre with glacial acetic acid and allowed to stand for 24 Hours before use. The acetic anhydride reacts with the water in perchloric acid and some traces in glacial acetic acid thereby making the resulting mixture practically anhydrous.

Standardisation of 0.1 N Perchloric acid

About 0.5g of Potassium hydrogen phthalate were weighed accurately and transfer in a 100 ml conical flask, previously dried at 120°C for 2 hours. 50ml of glacial acetic acid was added, warmed if necessary, to dissolve the salt completely. Cooled and titrated with 0.1M perchloric acid using crystal violet 0.1ml (2 drops) as an indicator. The end point is colour changes from violet to emerald green.

Each ml of 0.1M Perchloric acid = 0.02042g of Potassium hydrogen Phthalate

Assay of metronidazole:

About 3.0g of metronidazole sample was weighed accurately and transfer into a 250ml conical flask, 50ml of Glacial acetic acid was added, warmed gently. Cooled and titrated with 0.1M perchloric acid using α -Naphthol benzene as an indicator. The end point is colour changes from Blue to Blue-green.

Each ml of 0.1M perchloric acid = 0.02042g of Potassium hydrogen Phthalate.

Report:

The Molarity of 0.1M Perchloric acid=

The Percentage purity of given metronidazole was found to be=

Experiment No: 08

Assay of Dapsone


Aim:

To determine the percentage purity of given sample of tablet.

Principle:

Dapsone is a diamino diphenyl sulphone. It is assayed by direct diazotization titration. The free primary amino group present in dapsone is diazotized by nitrous acid and hydrochloric acid to form arene diazonium compound. The end point can be determined by using external indicator i.e. starch iodide paper, during titration with sodium nitrite an aromatic primary amine is diazotized. After the end point the drug solution containing a slight excess of nitrite comes in contact with starch iodide paper. It oxidizes iodide into iodine which gives blue colour with starch indicating the end point.

Reaction:



Procedure:

Preparation of 0.1M Sodium nitrite solution:

7.5g of sodium nitrite was dissolved in sufficient water to produced 1000ml.

Standardisation of 0.1M Sodium nitrite:

0.3g of sulphanilic acid was dissolved in 50ml of 2M hydrochloric acid, 3g of potassium bromide was added, cool in ice and titrate with 0.1M sodium nitrite solution using starch iodide paper as external indicator.

Each ml of 0.1M Sodium nitrite= 0.01732g of $C_6H_7NO_3S$.

Assay of Dapsone tablet:

20 tablets were weighed and powdered. Tablet powder equivalent to 0.25gm dapsone was weighed and dissolved in a mixture of 15ml of 2M hydrochloric acid. The solution was cooled to about 15°C and carry out sodium nitrite titration using starch iodide paper as an external indicator. End point is immediate appearance of blue colour.

Each ml of 0.1M sodium nitrite= 0.01241g of $C_{12}H_{12}N_2O_2S$.

Report:

The Molarity of 0.1M sodium nitrite=

The percentage purity of given dapsone tablet was found to be=

Experiment No:09

Assay of Chlorpheniramine Maleate



Aim:

To determine the percentage purity of given sample of Chlorpheniramine maleate.

Principle:

Chlorpheniramine maleate chemically known as 3-(4-chlorophenyl) N, N- Dimethyl-3-(2-pyridyl) propylamine hydrogen maleate is an anti-histamine H1 receptor antagonist. Mineral acid salts of weak nitrogen bases hydrolyse so extensively in aqueous or hydro-alcoholic solution that is possible to titrate the liberated acid with a strong mineral base. Titration of the maleate salt of the drug in water against the sodium hydroxide leads to the formation of water turbidity as the titration proceeds. To prevent this precipitation, alcohol has been used. Since alcohol is basic with respect to water as a solvent, dissolved bases react less strongly alkaline, their salts react more strongly acid, and the end points of the titrations are greatly sharpened, here an aqueous solution of Chlorpheniramine maleate was titrated with aqueous NaOH, the turbidity was formed with appearance of pink is the end point.

Reaction:



Procedure:**Preparation of 0.01M Sodium hydroxide:**

Weighed accurately about 0.4gm of Sodium hydroxide pellet in a clean 1000ml standard flask then completely dissolved with 100ml distilled water, and make up to 1000 ml with distilled water.

Standardisation of 0.01M sodium hydroxide:

Weighed accurately about 0.5gm of potassium hydrogen phthalate and transferred into 1000ml conical flask. Then add 75ml of distilled water to dissolve and titrated with 0.01M sodium hydroxide solution using 0.1ml phenolphthalein as an indicator. The end point is appearance of permanent pale pink colour.

Each ml of 0.01M Sodium hydroxide= 0.002042g of $C_8H_5KO_4$

Assay of Chlorpheniramine maleate

Twenty tablets were weighed and grind into a fine powder. An amount of powder equivalent to 200mg of chlorpheniramine maleate was weighed accurately into 100ml standard flask, 70ml of neutral alcohol was added and shaken for about 20 min. Then, the volume was made up to mark with neutral alcohol, mixed well and filtered using whatmann No 42 filter paper. The first 10 ml portion of the filtrate was discarded. An aliquot of the drug solution containing 2.0-20.0 mg of chlorpheniramine maleate was measured accurately and transferred into a clean 100 ml conical flask and the total volume was brought to 10 ml with neutral alcohol. Then 2 drops of 0.5% phenolphthalein indicator was added and the solution was titrated with standard 0.01M Sodium hydroxide solution until a permanent pink colour was obtained.

Each ml of 0.01M Sodium hydroxide= 0.002042g of $C_8H_5KO_4$.

Report:

The molarity of 0.01M Sodium hydroxide=



The percentage purity of given Chlorpheniramine maleate tablet was found to be=

Experiment: 10**Assay of Benzyl Penicillin****Aim:**


To determine the percentage purity of given sample of Benzyl penicillin tablet.

Principle:

Benzyl penicillin is assayed by iodometric titration method. The titration in which equivalent amount of iodine is liberated to form potassium iodide by the sample and the liberated iodine is titrated against standard sodium thiosulphate solution. This type of indirect titration determination of compounds is called iodometric titration. In this titration benzyl penicillin is first hydrolysed with sodium hydroxide solution converted to corresponding penicilloic acid (dicarboxylic acid). Then penicilloic acid further treated with mineral acid to form D-Penicillamine and benzyl penicillic acid. An obtained D-penicillamine is further oxidized quantitatively iodine to give disulphide, excess of iodine is back titrated with 0.02M sodium thiosulphate, equivalent amount of liberated iodine can be measured by titration with sodium thiosulphate using starch as an indicator, which is added near the end point as it get hydrolysed by hydrochloric acid and iodine gets trapped in the matrix of starch. Due to this there is no continuous liberation of iodine. An end point is blue to apple green.

Reaction:**Standardisation:**

Assay:



Procedure:

Preparation of 0.02N Sodium thiosulphate:

0.02N sodium thiosulphate (4.564 gm of sodium thiosulphate and 250 mg of sodium carbonate) dissolved in 100 ml of water then make up to 1000ml with water.

Standardization of 0.02N sodium thiosulphate:

Dissolved accurately weighed 0.2g of potassium bromate in 250ml of water taken in a conical flask. From this take 50 ml of the solution, add 2g of potassium iodide, 3ml of 2M HCl and titrate with sodium thiosulphate solution using starch as an indicator until the blue colour is discharged.

Each ml of 0.01N sodium thiosulphate= 0.002784g of KBr.

Assay of Benzyl penicillin:

Weighed accurately about 0.1gm of sodium salt of benzyl penicillin taken in a stoppered flask, dissolved in 10ml of water and dilute to 100ml. 10ml of solution was transferred into iodine flask, 5ml of 1N sodium hydroxide was added, allowed to stand for 20 minutes. Then freshly prepared buffer solution, 5ml of 1N hydrochloric acid and 25ml of excess of 0.02N iodine solution were added to the stoppered flask and kept aside for 20 minutes in a dark place. Excess of iodine is titrated with 0.02N sodium thiosulphate using freshly prepared starch solution as an indicator. The end point is discoloration of blue colour. To another 10ml of initial solution add 20ml of the buffer solution, allowed to stand for 20 minutes in the dark place and titrate with the same. The difference between two titration represents the volume of 0.02N iodine equivalent to the total amount of penicillin present in the given sample of benzyl penicillin.

Report:

The molarity of 0.01M sodium thiosulphate=

The percentage purity of given Benzyl penicillin tablet was found to be=

Experiment No: 11

Synthesis of Phenytoin from Benzil by Microwave Irradiation

Aim:


To synthesis and submit phenytoin from benzil by using microwave oven and calculate its percentage yield.

Principle:

The microwave region of Electro Magnetic spectrum light between IR irradiation and frequency corresponding to the wavelength of 1cm^{-1} to 1m. The organic compound can be heated by applying energy in the form of microwave high frequency IR irradiation. The use of this technique can have substantial saving time for laboratory synthesis of drugs and chemicals.

Phenytoin is prepared by condensation of benzil and urea under reflux condensation to give heterocyclic compound Pinnacol. When the pinnacol is treated with sodium hydroxide it's rearrangement to produce phenytoin as sodium salt and the esterification gives phenytoin as crude product. This reaction is completed to take place only 2 to 3 minutes microwave irradiation.

Reaction:



Chemicals Required:

- Benzil = 1.25g
- Urea = 0.75g
- NaOH (30%) = 3.75 ml
- Methanol = 18.75 ml

Procedure:

1.25gm of benzil, 0.75gm urea, 3.75ml 30% sodium hydroxide and 18.75 ml methanol was taken in a beaker and kept in a Microwave oven for 3 minutes. Then it was cooled to room temperature and 50ml of water was poured, mixed well and filtered. To the filtrate con. HCl, was added to get the crude product. Then the product was recrystallized from ethanol.

Melting Point: 295°C to 298°C

Use:

Anti-convulsant

Report:

Phenytoin was synthesised and reported the following

1. Theoretical yield =
2. Practical yield =
3. Percentage yield =

Experiment No: 12

Synthesis of Aspirin Assisted by Microwave Oven

Aim:


To Prepare and Submit aspirin from Salicylic acid and Calculate the Percentage Yield.

Principle:

The microwave region of Electro Magnetic spectrum light between IR irradiation and frequency corresponding to the wavelength of 1cm-1 to 1m. The organic compound can be heated by applying energy in the form of microwave high frequency IR irradiation. The use of this technique can have substantial saving time for laboratory synthesis of drugs and chemicals.

Aspirin is prepared from Salicylic acid, in this method salicylic acid react with acetic anhydride in the presence of sulphuric acid to form crystalline nature of aspirin.

Reaction:



Procedure:

Two reaction consist of a two-neck flask with magnetic stir bar, temperature, sensor and intensive cooler. A mixture of salicylic acid and acetic anhydride is filled in the reaction flask and three drop of conc. Sulphuric Acid are added the apparatus is installed by means of a glass tube in the microwave system. Standard refluxing apparatus for microwave system. The reaction mixture is heated under stirring for 90 sec. with 900 W to 140 C, during the following cooling down the clear yellowish solidified to a compact white crystalline mass

Melting Point: 135⁰C

Use:

Anti-inflammatory agent

Result:

Aspirin was synthesised and reported the following

1. Theoretical yield =
2. Practical yield =
3. Percentage yield =

Experiment No 13:

Drawing structure and Reaction using Chems sketch

CHEM SKETCH

ACD/ChemSketch is a molecular modeling program used to create and modify images of chemical structures. Also, there is a software that allows molecules and molecular models displayed in two and three dimensions, to understand the structure of chemical bonds and the nature of the functional groups.

The program offers some advanced features that allows the molecules rotate and apply colour to improve visualization. It has several templates with ions and functional groups with the possibility to add text and use other tools to optimize productions created by the software.

Application

Using ACD/ChemSketch is primarily for educational use. With this program it is possible to write and perform chemical equations, diagrams laboratories and chemical structures of various entity.

Computational chemistry software		
<u>Cheminformatics</u>	<u>Open-source</u>	<ul style="list-style-type: none"> ❖ <u>Avalon Cheminformatics Toolkit</u> ❖ <u>Bioclipse</u> ❖ <u>Blue Obelisk</u> ❖ <u>Chemistry Development Kit</u> ❖ <u>ECCE</u> ❖ <u>JOELib</u> ❖ <u>OELib</u> ❖ <u>Open Babel</u> ❖ <u>RDKit</u>
	<u>Proprietary</u>	<ul style="list-style-type: none"> ❖ <u>Canvas</u> ❖ <u>Chemicalize</u> ❖ <u>Discovery Studio</u>
<u>Chemical kinetics</u>	<u>Open-source</u>	<ul style="list-style-type: none"> ❖ <u>APBS</u> ❖ <u>Cantera</u> ❖ <u>KPP</u>
	<u>Proprietary</u>	<ul style="list-style-type: none"> ❖ <u>Autochem</u> ❖ <u>Chemical WorkBench</u> ❖ <u>CHEMKIN</u> ❖ <u>COSILAB</u> ❖ <u>DelPhi</u> ❖ <u>Khimera</u>

Molecular modelling visualization	<u>Open-source</u>	<ul style="list-style-type: none"> ❖ Ascalaph Designer ❖ Avogadro ❖ BALL ❖ Biskit ❖ CPMD ❖ Gabedit ❖ Ghemical ❖ Jmol ❖ Molekel ❖ PyMOL ❖ QuteMol ❖ RasMol
	<u>Proprietary</u>	<ul style="list-style-type: none"> ❖ Abalone ❖ ACD/ChemSketch ❖ Atomistix ToolKit ❖ ChemDraw ❖ ChemWindow ❖ EzMol ❖ Gaussian ❖ Maestro ❖ MarvinSketch ❖ MarvinView ❖ MODELLER ❖ MolecularOperating Environment ❖ SAMSON ❖ Spartan ❖ UCSF ❖ Chimera ❖ VMD
<u>Molecular docking</u>	<u>Open-source</u>	<ul style="list-style-type: none"> ❖ AutoDock ❖ AutoDock Vina ❖ FlexAID ❖ rDock
	<u>Proprietary</u>	<ul style="list-style-type: none"> ❖ Glide ❖ LeDock ❖ Molecular Environment <p align="right">Operating</p>
Molecular dynamics	<u>Open-source</u>	<ul style="list-style-type: none"> ❖ CP2K ❖ GROMACS ❖ LAMMPS ❖ OpenMM ❖ PLUMED
	<u>Proprietary</u>	<ul style="list-style-type: none"> ❖ Abalone ❖ AMBER ❖ CHARMM ❖ Desmond ❖ GROMOSNAMD

Quantum chemistry	<u>Open-source</u>	<ul style="list-style-type: none"> ❖ ABINIT ❖ ACES(CFOUR) ❖ AIMAll ❖ BigDFT ❖ CP2K ❖ DACAPO ❖ Dalton ❖ DP code ❖ FreeON ❖ HORTON ❖ MADNESS ❖ MPQC ❖ NWChem ❖ Octopus ❖ PARSEC ❖ PSI ❖ PyQuante ❖ PySCF ❖ Quantum ESPRESSO ❖ (PWscf) ❖ RMG ❖ SIESTA ❖ VB2000 ❖ YAMBO code
	<u>Proprietary</u>	<ul style="list-style-type: none"> ❖ ADF ❖ AMPAC ❖ DMol3 ❖ CADPAC ❖ CASINO ❖ CASTEP ❖ COLUMBUS ❖ CONQUEST ❖ CPMD ❖ CRUNCH ❖ CRYSTAL ❖ DIRAC ❖ Firefly ❖ GAMESS (UK) GAMESS (US) ❖ Gaussian ❖ Jaguar ❖ MOLCAS ❖ MOLPRO ❖ MOPAC ❖ ONETEP ❖ OpenAtom ❖ ORCA ❖ PLATO ❖ PQS

		<ul style="list-style-type: none"> ❖ Q-Chem ❖ Quantemol ❖ Scigress ❖ Spartan ❖ TeraChem ❖ TURBOMOLE ❖ VASP ❖ WIEN2k ❖ XMVB
Skeletal Structure drawing	<u>Open-source</u>	<ul style="list-style-type: none"> ❖ BKChem ❖ JChemPaint ❖ JME Molecule Editor ❖ Molsketch ❖ XDrawChem
	<u>Proprietary</u>	<ul style="list-style-type: none"> ❖ ACD/ChemSketch ❖ ChemDoodle ❖ ChemDraw ❖ ChemWindow ❖ MarvinSketch
Others		<ul style="list-style-type: none"> ❖ Aqion ❖ Eulim ❖ EXC code ❖ GenX ❖ GSim ❖ ICM (ICM-Browser) ❖ Materials Studio ❖ Molden ❖ OpenChrom ❖ RubyChem ❖ SASHIMI